


## New 9600-10D PUREIL pH electrode suitable for accurate and fast measurement of low conductivity water

Innovative pH electrode equipped with newly developed gelled ionic liquid for the liquid junction of the reference electrode. Unlike conventional pH electrodes, it does not leak potassium chloride (KCl), which can affect the measurement or pH value of low ionic strength aqueous samples and certain critical solutions making it an excellent choice for such applications.


- ✓ Ideal for low conductivity water ( $< 100 \mu\text{S}/\text{cm}$ ) such as pure water samples (e.g. distilled water, deionized water, rain water, reverse osmosis water, boiler water)
- ✓ Perfect for solutions that react with KCl or contain ions that form precipitates with  $\text{K}^+$  or  $\text{Cl}^-$  such as drug injection solution
- ✓ Only slight dissolution of ionic liquid in sample during measurement
- ✓ Maintenance-free and waterproof
- ✓ Integrated with temperature sensor for automatic temperature compensation
- ✓ Comes with sensor holder for mounting to electrode stand

### 9600-10D

### PUREIL



Sensor holder



Conventional pH Electrode

9600-10D PUREIL pH Electrode

3.33M KCl diffuses into the sample

Only slight dissolution of ionic liquid [ $< 1/10000$  of conventional pH electrode in terms of concentration]

Traditional Salt Bridge

Ionic Liquid Salt Bridge



This product is the result of joint development with Kyoto University as part of the Japan Science and Technology Agency (JST) advanced measurement analysis technology and equipment development program.

Conventional pH electrode utilizes reference electrode with KCl solution as salt bridge electrolyte. KCl leaks into the sample through the liquid junction to complete the electrical circuit during pH measurement. Diffusion of KCl dominates the liquid junction potential.

» In testing low ionic strength aqueous solutions, KCl can contaminate sample resulting in increase in conductivity and change in pH value. Thus, it is difficult to achieve accurate pH measurement in these solutions using conventional pH electrode with KCl salt bridge.

Modern pH electrode like the 9600-10D PUREIL utilizes reference electrode with new salt bridge called ionic liquid (room-temperature molten salt) which serves as point of contact between the reference electrolyte and sample. The cations (C<sup>+</sup>) and the anions (A<sup>-</sup>) of the ionic liquid partitioned in the sample determine the liquid junction potential difference between reference electrolyte and sample.


» Sample is not exposed to KCl. Therefore, no KCl contamination. The 9600-10D PUREIL pH electrode shows very stable liquid junction potential, faster response, and accurate pH values in low ionic strength aqueous solutions.

## Measurement Examples

### Low Conductivity Water

Since there is no concentrated KCl diffusion from the pH electrode, the sample pH value remains unaffected during measurement. Also, since the ionic liquid dissolves quickly, the liquid junction potential stabilizes fast leading to a more reliable and stable pH reading.

#### Dilute Sulfuric Acid: 20 μmol/dm<sup>3</sup> (100 μS/cm)




### Drug Injection Solution

Solutions that form precipitates with KCl such as drug injection solution can be measured without any adverse reaction.

#### Conventional pH electrode

Precipitate forms when sample reacts with KCl.



#### 9600-10D

No formation of precipitate.

Please note that the 9600-10D PUREIL pH electrode is not suitable for measuring samples that contain surfactant, heavy metal ions, and organic solvents. For acidic samples, citrate buffer (pH 3.776 at 25°C) is recommended for calibration instead of phthalate buffer (pH 4.01 at 25°C) to achieve a more accurate measurement [ $< 0.034$  pH error].

## PUREIL Electrode Specifications



<b>Model</b>	9600-10D
<b>pH Measurement Range</b>	0 – 14
<b>Operating Temperature Range</b>	0 – 40 °C
<b>Storage Temperature Range</b>	0 – 50 °C
<b>Internal Electrode</b>	Ag/AgCl
<b>Liquid Junction</b>	PVDF-HFP, ionic liquid
<b>Liquid Junction Height</b>	Approx. 16mm from the electrode tip
<b>Wetted Materials</b>	Glass, PBT, PVDF-HFP, ionic liquid
<b>Electrode Dimensions</b>	210 mm (L) x 12 mm (Ø)
<b>Cable</b>	1 m

**HORIBA Instruments (Singapore) Pte. Ltd.**  
83 Science Park Drive, #02-02A, The Curie, Singapore 118258  
Phone: 65 6908-9660  
Fax: 65 6745-8155  
[www.horiba-laqua.com](http://www.horiba-laqua.com)  
e-mail: [laqua@horiba.com](mailto:laqua@horiba.com)

**HORIBA UK Limited**  
Kyoto Close, Moulton Park, Northampton NN3 6FL  
Phone: 44 (0) 1604 542567  
Fax: 44 (0) 1604 542699  
[www.horiba.com/uk](http://www.horiba.com/uk)  
e-mail: [waterquality@horiba.com](mailto:waterquality@horiba.com)

**HORIBA Instruments Incorporated**  
9755 Research Drive, Irvine, California 92618 USA  
Phone: +1 949 250 4811  
FAX: +1 949 250 0924  
[www.horiba.com/us/en](http://www.horiba.com/us/en)  
e-mail: [labinfo@horiba.com](mailto:labinfo@horiba.com)



[www.horiba-laqua.com](http://www.horiba-laqua.com)

**IMS**  
HORIBA Group is operating Integrated Management System (IMS)  
ISO9001 JOA-0298 / ISO14001  
JOA-E-90039 / ISO13485  
JOA-MD0010 / OHSAS18001 JOA-OH0068Rorum perid novis esimoente,

