

TEST REPORT
IEC 61010-1
Safety requirements for electrical equipment for measurement,
control, and laboratory use
Part 1: General requirements

Report Number: 64.111.21.00851.01 Rev.00

Date of issue: 2021-09-13

Total number of pages.....: 87

Name of Testing Laboratory: TÜV SÜD Certification and Testing (China) Co., Ltd. Guangzhou
preparing the Report.....: Branch

Applicant's name.....: Qingdao Haier Biomedical Co., Ltd.

Address: Haier Industrial Park, Economic Technology Development Zone,
Qingdao, 266510, P.R. China

Test specification:

Standard: IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016

Test procedure: CE LVD

Non-standard test method: N/A

TRF template used: IECEE OD-2020-F1:2020, Ed.1.3

Test Report Form No......: IEC61010_1P

Test Report Form(s) Originator: VDE Prüf- und Zertifizierungsinstitut GmbH

Master TRF: 2021-04-12

Copyright © 2021 IEC System of Conformity Assessment Schemes for Electrotechnical Equipment and Components (IECEE System). All rights reserved.


This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.

This report is not valid as a CB Test Report unless signed by an approved IECEE Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.

General disclaimer:

The test results presented in this report relate only to the object tested.

This report shall not be reproduced, except in full, without the written approval of the Issuing NCB. The authenticity of this Test Report and its contents can be verified by contacting the NCB, responsible for this Test Report.

Test item description :	Drying Oven	
Trade Mark :	Haier	
Manufacturer	Same as applicant	
Model/Type reference :	HZS-60, HFS-160	
Ratings :	220-240V~, 50/60Hz, class I, IPX0, 900W for HZS-60, 2500W for HFS-160	
Responsible Testing Laboratory (as applicable), testing procedure and testing location(s):		
<input checked="" type="checkbox"/>	Testing Laboratory:	TÜV SÜD Certification and Testing (China) Co., Ltd. Guangzhou Branch
	Testing location/ address	Location 1: 5F, Communication Building, 163 Pingyun Rd, Huangpu Ave. West Guangzhou 510656 P. R. China Location 2: B1 Floor, No. 3 Chuangqi Building, TusPark, 63 Chuangqi Road, Shilou Town, Panyu District, Guangzhou 511447, P.R. China
	Tested by (name, function, signature) :	Niu Weicong, project handler <i>Niu Weicong</i>
	Approved by (name, function, signature) .. :	Wang Xujun, designated reviewer <i>Wang Xujun</i>
		
<input type="checkbox"/>	Testing procedure: CTF Stage 1:	N/A
	Testing location/ address	N/A
	Tested by (name, function, signature) :	N/A
	Approved by (name, function, signature) .. :	N/A
<input type="checkbox"/>	Testing procedure: CTF Stage 2:	N/A
	Testing location/ address	N/A
	Tested by (name + signature) :	N/A
	Witnessed by (name, function, signature) .. :	N/A
	Approved by (name, function, signature) .. :	N/A
<input type="checkbox"/>	Testing procedure: CTF Stage 3:	N/A
<input type="checkbox"/>	Testing procedure: CTF Stage 4:	N/A
	Testing location/ address	N/A
	Tested by (name, function, signature) :	N/A
	Witnessed by (name, function, signature) .. :	N/A
	Approved by (name, function, signature) .. :	N/A
	Supervised by (name, function, signature) :	N/A

List of Attachments (including a total number of pages in each attachment)		
Document No.	Documents included / attached to this report (description)	Page No.

